

Larident Srl

Via Lamaneigra, 12 r (Ponte Rimassa) - 16030 Tribogna Genova

Tel. +39 0185 939 481 - Fax +39 0185 938 595

info@larident.it - www.larident.it



TECHNICAL SHEET

PRODUCT NAME	ECONOMICAL FACE SHIELD WITH 6 VISORS.
CODE	R5
COMPONENT	PROTECTIVE VISOR: POLICARBONATO ANTIAPPANNANTE. VISOR THICKNESS: 0.178
THERMAL CHARACTERISTIC	THE VISOR IS REUSABLE AND NOT STERILISABLE
PHYSICAL CHARACTERISTICS, MECHANICAL CHARACTERISTICS AND PACKING	PACKING IN CARTON BOX: 6 PCS THE VISOR COVER NOSE MOUTH AND EYES THE PROTECTIVE VISOR HAS A DOUBLE PROTECTIVE COATING THE VISOR IS VERY TRANSPARENT
CHEMICAL CHARACTERISTIC	IT RESISTS TO THE CHEMICAL DISINFECTANTS.
SAFETY DATA SHEET	NOT APPLICABLE
PRECAUTION	WARNING : USE ONLY AS FACE SHIELD AGAINST SPLASHING LIQUID OF WATER SALIVA.
INSTRUCTIONS FOR USE	SUPPLIED IN "NON-STERILE" CONDITION. AUTOCLAVABLE. THE PRODUCT IS SINGLE USE AVAILABLE INSTRUCTION FOR FITTING INSIDE CARTON BOX.
ANALYSIS OF THE RISKS REF.D.L.46 24/2/97	AVAILABLE ONLY IN ITALIAN
CLAS. REGULATION 93/42 CEE	EN166, cat. II