

3M Science.
Applied to Life.™

3M™ Scotchbond™ Universal Plus
Adhesive

**More than just
a strong bond.**





3M™
Scotchbond™
Universal Plus
Adhesive takes
universal adhesives to
the next level.


It keeps all the benefits of the original 3M™ Scotchbond™ Universal Adhesive, the first truly universal adhesive, while adding even greater control and predictability – because a good adhesive should provide more than just a strong bond.

From excellent...

Built on the clinically proven and world's most researched universal adhesive.

 **All direct and indirect indications**
light-curing, self-curing or dual-curing
methacrylate-based composites / cements

 **All etching techniques**
total-etch, self-etch, and selective
enamel etch

 **Bonds to all restorative materials**
without the need for separate primers

 **Available in vial and unit dose**
for efficient hygiene management

 **Virtually no post-operative sensitivity**
in field evaluation:

99% of dentists reported
that 3M™ Scotchbond™ Universal Plus
Adhesive shows virtually no post-
operative sensitivity.*

...to exceptional

Formulation carefully optimised to achieve unique features.

 **First radiopaque universal adhesive**
reducing the risk of X-ray misdiagnosis and
invasive overtreatment*


 **Bonds and seals caries-affected dentine**
maximising preservation of natural tooth structure,
supporting minimally invasive preparation
guidelines*

 **Gold standard adhesion to all dental substrates**
including glass ceramics

 **Full dual- and self-cure compatibility**
no need for an additional dual-cure activator bottle

 **BPA free formulation**
free of BPA derivatives like
BisGMA




 **Fully aligned system with 3M™ RelyX™
Universal Resin Cement**
excellent bond strength for virtually all resin
cement indications.* Adhesive cured by cement.

*Field evaluation with 309 dentists and over 20,000 applications



Added radiopacity


On average, dentists see at least one hard to interpret X-ray per day due to radiolucency under an existing restoration. This could lead to misdiagnosis and invasive overtreatment. 3M™ Scotchbond™ Universal Plus Adhesive is the first radiopaque universal adhesive with dentine-like radiopacity, helping to overcome this issue.*



Non-Radiopaque Adhesive

3M™ Scotchbond™ Universal Plus Adhesive

Intentional adhesive pooling under Class I restorations – X-rays of extracted human molars, filling material 3M™ Filtek™ One Bulk Fill Restorative. 3M internal data, H. Loll et al., J Dent Res 99 (Spec Iss A): No. 757, 2020



“Radiolucent adhesive pooling underneath a restoration doubles the risk of restoration replacement.”

Dr. John O. Burgess,
University of Alabama, USA

Ordering information

Item no. Product information

3M™ Scotchbond™ Universal Plus Adhesive – Vial

41293	Intro Kit Vial 1 Vial 5ml, 50 Disposable Brushes, 1 Etchant Syringe 3ml, 25 Dispensing Tips
41294	Refill Vial – 1 × 5ml
41295	Refill Vial – 3 × 5ml
41263	3M™ Scotchbond™ Universal Etchant Etching Gel 2 Syringes x 3ml, 50 Dispensing Tips



Vial
with smooth and even surface –
easy to wipe off excess

3M™ Scotchbond™ Universal Plus Adhesive – Unit Dose

41297	Intro Kit Unit Dose 50 unit doses, 1 Etchant Syringe 3ml, 25 Dispensing Tips
41298	Unit Dose 100 100 unit doses, 1 Etchant Syringe 3ml, 25 Dispensing Tips
41299	Unit Dose 200 200 unit doses, 1 Etchant Syringe 3ml, 25 Dispensing Tips
41304	Unit Dose 400 400 unit doses, 1 Etchant Syringe 3ml, 25 Dispensing Tips



Unit Dose
for efficient hygiene management –
reduces the risk of cross
contamination

*3M internal data

3M Oral Care
3M United Kingdom PLC
Charnwood Campus
10 Bakewell Road
Loughborough
Leicestershire, LE11 5RB
+0800 626 578

3M Ireland Ltd
The Iveagh Building
Carrickmines Park
Carrickmines Dublin 18
Ireland
+353 1 280 3555



[3M.co.uk/Scotchbond-plus](https://www.3m.co.uk/Scotchbond-plus)

☎ 08705 360036

✉ OralCareUK@mmm.com

📱 @3MOralCareUK

3M, Filtek, RelyX and Scotchbond are trademarks of 3M Company. © 3M 2020. All rights reserved. G0113A / 17422317_UK20201111